



Medion Diagnostics

New

.... DAC III

**Cellwashing centrifuge for small and
midsize labs**



➤ **characteristics**

- benchtop compact centrifuge
- individual programming for centrifugation time, wash cycles, and NaCl volume
- convenient touch pad programming
- ultra-quiet brushless drive

➤ **capacity**

- 12-place removable rotor
- or
- 24-place removable rotor

➤ **price**

- upon request

Further information available at:

Medion Diagnostics GmbH, Bonnstrasse 9, CH-3186 Dürdingen

Phone: +41 (26) 492 85 11; Fax: +41 (26) 492 86 00

Customer Service: Phone: +41 (26) 492 87 02; Fax +41 (26) 492 86 56